



Solvent Exchange from Dichloromethane to n-Hexanes Using the DryVap® Concentrator System

Stephen J MacDonald, Horizon Technology, Ins., Salem, NH

Introduction

The analysis of pesticides for environmental extracts typically requires a solvent exchange from dichloromethane (DCM) to an ECD-compatible solvent such as n-hexane. The prescribed method of solvent exchange and concentration is given in EPA Method 3510 and uses a Kuderna-Danish (K-D) apparatus achieve the results.

The K-D method is very time consuming and can take up to 1.5 hr to complete. It is also sensitive to human error and mistakes such as evaporating too far or too fast or not conditioning the sodium sulfate correctly can often result in lost samples.

The DryVap® concentrator system can automatically dry and concentrate the extract in as little as 30 min. Solvent exchange can be performed by adding 30 mL hexane to the concentrate and repeating the automated process.

The DryVap® concentrator system eliminates all but one of the manual steps in the process and accomplishes the task faster than manual K-D evaporations. Furthermore, the final dried, solvent-exchanged, and concentrated sample is kept at room temperature under a blanket of inert gas and can be allowed to stand upon completion on the DryVap®.

To demonstrate the performance of the DryVap® system, three samples were spiked with a pesticides mixture and were then dried, concentrated, and exchange.

Instrumentation

- Horizon Technology
 - DryVap® Concentrator System
 - DryDisk®
- Agilent 5890 GC-ECD

Method Summary

1. Spike three 200 mL aliquots of DCM with 20 ng of organochlorine pesticides mix.
2. Pour the samples into the DryDisk® station of the DryVap® System.
3. Start the DryVap®, the samples will be dried and concentrated automatically to 0.9 mL.
 - Main Vacuum: 10 in. Hg
 - Nitrogen Pressure: 20 psi
 - Dry Volume: 200 mL
 - Heater Power: 5
 - Auto Rinse: 0
 - Spurge Heater: OFF
4. Add 30 mL of n-hexane to the concentrator tube.
5. Using a Pasteur pipette, mix the concentrated extract in the nipple with the exchange solvent.

Results

Table 1 shows the average recoveries and the relative standard deviation (RSD) for the three completed runs. The automation of the DryVap® system allows for the RSD to be less than 5 % (with the exception of three compounds). Also, the high recovery evident for both 4,4'-DDT and Endrin indicate that there is no degradation of thermally labile compounds.

Table 1: Average Recoveries and RSD for Three Organochlorine Pesticide Spikes

Compound	Average Rec.	RSD
a-BHC	78	4.8
Lindane	81	4.7
b-BHC	85	5.3
d-BHC	85	4.2
Heptachlor	87	4.3
Aldrin	79	3.6
Heptachlor epoxide	80	3.3
g-Chlordane	87	9.3
a-Chlordane	79	3.2
4,4'-DDE	85	0.7
Endosulfan I	77	2.7
Dieldrin	78	2.0
Endrin	103	2.9
4,4'-DDD	85	1.8
Endosulfan II	82	1.9
4,4'-DDT	92	1.3
Endrin aldehyde	71	3.7
methoxychlor	89	3.0
endosulfan sulfate	82	3.7
Endrin ketone	79	3.2

Conclusions

The DryVap[®] concentration system allows for controlled drying, concentration, and solvent exchange to be completed in very little time. It removes the human element from the process more than conventional methods. This allows for fewer mistakes to be made by the analyst and a higher rate of reproducibility, both of which can save a company time and money.

EUROPE

Main Office: +46 18 565900
Toll Free: +800 18 565710
Fax: +46 18 591922
Order Tel: +46 18 565710
Order Fax: +46 18 565705
order@biotage.com
Support Tel: +46 18 56 59 11
Support Fax: +46 18 56 57 11
eu-1-pointsupport@biotage.com

NORTH & LATIN AMERICA

Main Office: +1 704 654 4900
Toll Free: +1 800 446 4752
Fax: +1 704 654 4917
Order Tel: +1 704 654 4900
Order Fax: +1 434 296 8217
ordermailbox@biotage.com
Support Tel: +1 800 446 4752
Outside US: +1 704 654 4900
us-1-pointsupport@biotage.com

HORIZON TECHNOLOGY

Tel: +1 603 893 3663
Toll Free: +1 800 997 2997
Fax: +1 603 893 4994
sales@biotage.com

JAPAN

Tel: +81 3 5627 3123
Fax: +81 3 5627 3121
jp_order@biotage.com
jp-1-pointsupport@biotage.com

CHINA

Tel: +86 21 68162810
Fax: +86 21 68162829
cn_order@biotage.com
cn-1-pointsupport@biotage.com

KOREA

Tel: +82 31 706 8500
Fax: +82 31 706 8510
korea_info@biotage.com
kr-1-pointsupport@biotage.com

INDIA

Tel: +91 22 4005 3712
india@biotage.com

Distributors in other regions are listed on www.biotage.com

Part Number: AN048-HOR

© 2009 Biotage. All rights reserved. No material may be reproduced or published without the written permission of Biotage. Information in this document is subject to change without notice and does not represent any commitment from Biotage. E&OE. A list of all trademarks owned by Biotage AB is available at www.biotage.com/legal. Other product and company names mentioned herein may be trademarks or registered trademarks and/or service marks of their respective owners, and are used only for explanation and to the owners' benefit, without intent to infringe.



a part of

