

Extraction of Barbiturates from Oral Fluid Using ISOLUTE® SLE+ After Collection with the NeoSal™ Collection Device Prior to GC/MS Analysis




Figure 1. Structure of Butabarbital.

Introduction

This application note describes the extraction of barbiturates from oral fluid matrix collected using the NeoSal collection devices, prior to GC/MS analysis

ISOLUTE® SLE+ supported liquid extraction plates and columns offer an efficient alternative to traditional liquid-liquid extraction (LLE) for bioanalytical sample preparation, providing high analyte recoveries, no emulsion formation, and significantly reduced sample preparation.

This application note describes an effective and efficient ISOLUTE SLE+ protocol optimized for 1 mL sample volumes. The simple sample preparation procedure delivers clean extracts and analyte recoveries greater than 88% with RSDs lower than 6% for all analytes.

Analytes

Butalbarbital, Butabarbital, Amobarbital, Pentobarbital, Secobarbital, Hexobarbital and Phenobarbital.

Sample Preparation Procedure

Format:

ISOLUTE SLE+ 1 mL Sample Volume Column, Part Number 820-0140-C

Sample Pre-treatment:

Following saliva collection, add 18 µL concentrated ammonium hydroxide to the matrix buffer contained in each collection device.

Sample Loading:

Load 1 mL of the pre-treated oral fluid device buffer matrix onto the column and apply a pulse of vacuum or positive pressure (3–5 seconds) to initiate flow. Allow the sample to absorb for 5 minutes.

Analyte Extraction:

Apply 2.5 mL of DCM/IPA (95/5, v/v) and allow to flow under gravity for 5 minutes. Apply a further aliquot of DCM/IPA (95/5, v/v, 2.5 mL) and allow to flow for another 5 minutes under gravity. Apply vacuum or positive pressure to pull through any remaining extraction solvent (5–10 seconds).

Post Elution and Reconstitution:

Dry the extract in a stream of air or nitrogen using a TurboVap® (10 psi at 40 °C for 40 mins).

Reconstitute the extracts with 250 µL ethyl acetate and vortex for 20 seconds before transferring to high recovery GC vials.

Dry the extract in a stream of air or nitrogen using a SPE Dry (40 °C, 20 to 40 L/min) or TurboVap® (10 psi at 40 °C for 40 mins).

Upon dryness, reconstitute with 80 µL ethyl acetate and 20 µL TMAH (trimethylanilinium hydroxide, 0.2M) and vortex for 20 seconds.

GC Conditions

Instrument

Agilent 7890A with QuickSwap

Column

Agilent J&W DB-5ms, 30 m x 0.25 mm ID x 0.25 µm

Carrier

Helium 1.2 mL/min (constant flow)

Inlet

260 °C, Splitless, purge flow: 50 mL/min at 1.0 min

Injection

1 µL

Wash Solvents

Methanol and ethyl acetate

Oven

Initial temperature 120 °C, hold for 1 minute

Ramp 12 °C/min to 192 °C

Ramp 100 °C/min to 330 °C, hold for 0.5 minutes

Post Run

Backflush for 1.6 minutes (2 void volumes)

Transfer Line

280 °C

MS Conditions

Instrument

Agilent 5975C

Source

230 °C

Quadrupole

150 °C


MSD mode

SIM

SIM Parameters

Table 1. Ions acquired in the Selected Ion Monitoring (SIM) mode.

SIM Group	Analyte	Target (Quant) Ion	1st Qual Ion	2nd Qual Ion
1	Butalbarbital	196	195	181
1	Butabarbital	169	184	211
2	Amobarbital	169	184	225
3	Pentobarbital	169	184	225
4	Secobarbital	196	195	181
5	Hexobarbital	235	81	169
6	Phenobarbital	232	146	175



Results

This optimized ISOLUTE® SLE+ protocol demonstrated analyte recoveries ranging from 102–107% as shown in **Figure 2**. RSDs were below 6% for all analytes.




Figure 2. Typical analyte % extraction recoveries (n=7) using the ISOLUTE® SLE+ protocol.




Figure 3. Typical ISOLUTE® SLE+ procedure.




Figure 4. Calibration curves for extracted levels of spiked oral fluid after collection with NeoSal devices, using 1 mL ISOLUTE® SLE+ format. Analyte concentrations spiked into each device are 5, 10, 25, 100 and 250 ng/mL showing r^2 values of 0.995 to 0.999.




Figure 5. GC/MS chromatography for NeoSal collected oral fluid spiked at 250 ng/mL. Early eluting peaks are visually poor in shape due to the acquisition of 6 SIMs but mass spectrometry is able to determine the quantification m/z with no contribution or interference.

Reagent Preparation

Concentrated ammonium hydroxide: Concentrated stock used to modify pH prior to extraction is commercially available (28–30%).

Table 2. Lower Limits of Quantitation (LLOQ) using NeoSal™ devices prior to optimized ISOLUTE® SLE+ procedure

Drug Analyte	LLOQ (ng/mL)
Butalbarbital	5
Butabarbital	5
Amobarbital	5
Pentobarbital	5
Secobarbital	5
Hexobarbital	5
Phenobarbital	5

Ordering Information

Part Number	Description	Quantity
820-0140-C	ISOLUTE® SLE+ 1 mL Sample Volume Column	30
PPM-48	Biotage® PRESSURE+ 48 Positive Pressure Manifold	1
SD-9600-DHS-EU	Biotage® SPE Dry 96 Sample Concentrator System 220/240V	1
SD-9600-DHS-NA	Biotage® SPE Dry 96 Sample Concentrator System 100/120V	1
C103199	TurboVap® LV Evaporator	1

EUROPE

Main Office: +46 18 565900
Toll Free: +800 18 565710
Fax: +46 18 591922
Order Tel: +46 18 565710
Order Fax: +46 18 565705
order@biotage.com
Support Tel: +46 18 56 59 11
Support Fax: +46 18 56 57 11
eu-1-pointsupport@biotage.com

NORTH & LATIN AMERICA

Main Office: +1 704 654 4900
Toll Free: +1 800 446 4752
Fax: +1 704 654 4917
Order Tel: +1 704 654 4900
Order Fax: +1 434 296 8217
ordermailbox@biotage.com
Support Tel: +1 800 446 4752
Outside US: +1 704 654 4900
us-1-pointsupport@biotage.com

JAPAN

Tel: +81 3 5627 3123
Fax: +81 3 5627 3121
jp_order@biotage.com
jp-1-pointsupport@biotage.com

CHINA

Tel: +86 21 2898 6655
Fax: +86 21 2898 6153
cn_order@biotage.com
cn-1-pointsupport@biotage.com

To locate a distributor,
please visit our website
www.biotage.com

Part Number: AN868

© 2016 Biotage. All rights reserved. No material may be reproduced or published without the written permission of Biotage. Information in this document is subject to change without notice and does not represent any commitment from Biotage. E&OE. A list of all trademarks owned by Biotage AB is available at www.biotage.com/legal. Other product and company names mentioned herein may be trademarks or registered trademarks and/or service marks of their respective owners, and are used only for explanation and to the owners' benefit, without intent to infringe.