

Cannabinol (CBN) Purification Using Flash Chromatography

Summary

Most of the early interest in medicinal cannabis applications focused on cannabidiol (CBD), and this is now especially true with various introductions and revisions to global regulatory and legislative standards on this topic. For example, in the US, with the passage of the 2018 Farm bill, some cannabis products were de-scheduled from the Controlled Substances Act. There has also been concurrent and increasing scientific interest in medical applications involving the so-called minor cannabinoids such as CBN, CBG, and CBC. This application note focusses on CBN.

Introduction

CBN is a degradation product of THC and is expressed naturally only in small quantities in cannabis plants. Through oxidative processes (light/heat/air¹ or bio-enzymatically²), THC oxidizes to CBN. Conversion of CBD to CBN is also possible using iodine (Federica Pollastro, 2018).

Whether formed by natural degradation/oxidation or by chemical/enzymatic reaction, there is always contamination with other cannabinoids, reagents, and by-products all which need removal in order to isolate a pure product.

Methods - Synthetic CBN

Flash system

Biotage® Isolera Dalton 2000 (APCI)

Column

Biotage® Sfär C18, 12-gram

Solvent

A: Water

B: Methanol

Equilibration

80%B for 2 CV

Gradient

80%B for 6 CV

100%B for 6 CV

Flow Rate

25 mL/min

Detection

UV @ 283 nm, 218 nm (collect)

UV 200–400 nm (monitor)

Mass monitor m/z -309, -313, +311, +315

Evaporation

Biotage® V-10 Touch

Methods - Natural CBN

Biotage® Isolera Dalton 2000 (APCI)

Column

Biotage® Sfär C18, 12-gram

Solvent

A: Water

B: Methanol

Equilibration

80%B for 2 CV

Gradient

80%B for 20 CV

Flow Rate

25 mL/min

Detection

UV @ 283 nm, 254 nm (monitor)

UV 200–350 nm (collect)

Mass monitor m/z +311, +315, +287, +317

Protocol

Synthetic CBN

- » Transfer 5 mL of crude CBN reaction mixture into a tared 20-mL scintillation vial
- » Evaporate solvent using a hi-boil method with V-10 Touch
- » Add methanol, QS to 5 mL (concentration 670 mg/mL)
- » Create flash chromatography method
- » Operate method
 - » Load 1 mL (670 mg) onto column when prompted

Natural CBN

- » To 100 mg of a 5-year old cannabis extract add 1 mL of methanol and dissolve completely
- » Using the same column as above, inject 0.1 mL

Results and Discussion

Synthetic CBN

The ability to scale purification is important in a production environment so chromatographic methods need to be simple, robust and efficient. The method utilized in this application is a simple step gradient that separates CBN from its reaction mixture impurities, Figure 1.




Figure 1. Simple step gradient purification method for synthetic CBN using both UV and mass detection. CBN is the big peak in the middle with a detected negative ionization mass of -309.




Figure 2. UV spectral analysis of purified CBN at the beginning of the peak (left) and at the end of the peak (right) shows the collected CBN fractions are pure. The amplitudes differ based on CBN concentration at the measurement point.

Both UV and mass detection were employed but only the UV signal was used to trigger fractionation (283 nm) leading to method simplicity. Mass detection was used to identify CBN as it has a unique molecular mass (310), which ionizes in the mass detector by adding a hydrogen (positive ionization) to create a detectable mass of +311 and subtracting a hydrogen (negative ionization) creating a detectable mass of -309. Most other cannabinoids have a molecular weight of 314.

Another useful technique to identify CBN is photo diode array UV (PDA), available on Biotage Isolera and Selekt systems. With PDA all the wavelengths emitted by the UV lamp are focused through the column effluent and the amount of light absorption at each wavelength is recorded and displayed in real time. The data can be analyzed post-run to evaluate UV peak purity and identify an individual compound, especially powerful when the compound's spectral footprint is known.

In the case of CBN, a literature search uncovered a paper showing both UV and mass spectra for various cannabinoids (Hazekamp, 2005). In this article, the CBN UV absorption spectra shows strong UV maxima at both 220 nm and 285 nm. By comparison, most other cannabinoids have strong UV maxima around 205 nm, medium absorbance near 230 nm, and minor absorption in the 270-280 nm range³.

For the synthetic CBN purification, the UV spectra shows the same UV maxima as shown in the Hazekamp article. Because of the high CBN concentration, the UV spectra was evaluated at the beginning and end of the peak, Figure 2. Because the UV maxima are the same at the front and rear, this suggests the collected fractions are pure CBN.

Post purification, the CBN fractions were evaporated using the V-10 Touch and a clear, semi crystalline product obtained, Figure 3.



Figure 3. Purified CBN from a synthetic reaction mixture.

Natural CBN

Because CBN is an oxidative degradation product of THC, its purification is a bit more challenging as the two cannabinoids elute closely to each other. To best facilitate this separation a simple isocratic method was employed using UV for fractionation and mass for identification. This method elutes CBN in 6 minutes with a detectable mass of +311; THC closely follows at ~7 minutes with a detectable mass of +315, Figure 4.




Figure 4. Isocratic purification of natural CBN created through THC oxidation. CBN is the large peak at 5 minutes with a detected mass of +311. THC elutes just after CBN and has a detected mass of +315.

As with the synthetic CBN, the UV spectra is a confirmatory tool that verified the compound identification, Figure 5.

Conclusion

Reversed-phase flash chromatography with diode-array UV and/or mass detection has been shown to both purify and identify CBN from synthetic and natural sources.

References

1. Reddit.com
2. US Patent US10538790B2
3. Hazekamp, A; Peltenburg, A; Verpoorte, R. Journal of Liquid Chromatography & Related Technologies®, 2005, 28, 2361–2382.




Figure 5. Purified degraded THC extract UV spectra identifies the tall green peak as CBN with UV maxima at 218 nm and 283 nm (left) and the adjacent peak to the right as THC with UV maxima of 209 nm and 278 nm (right).

EUROPE

Main Office: +46 18 565900
 Toll Free: +800 18 565710
 Fax: +46 18 591922
 Order Tel: +46 18 565710
 Order Fax: +46 18 565705
 order@biotage.com
 Support Tel: +46 18 56 59 11
 Support Fax: + 46 18 56 57 11
 eu-1-pointsupport@biotage.com

NORTH & LATIN AMERICA

Main Office: +1 704 654 4900
 Toll Free: +1 800 446 4752
 Fax: +1 704 654 4917
 Order Tel: +1 704 654 4900
 Order Fax: +1 434 296 8217
 ordermailbox@biotage.com
 Support Tel: +1 800 446 4752
 Outside US: +1 704 654 4900
 us-1-pointsupport@biotage.com

JAPAN

Tel: +81 3 5627 3123
 Fax: +81 3 5627 3121
 jp_order@biotage.com
 jp-1-pointsupport@biotage.com

CHINA

Tel: +86 21 68162810
 Fax: +86 21 68162829
 cn_order@biotage.com
 cn-1-pointsupport@biotage.com

KOREA

Tel: + 82 31 706 8500
 Fax: + 82 31 706 8510
 korea_info@biotage.com
 kr-1-pointsupport@biotage.com

INDIA

Tel: +91 22 4005 3712
 india@biotage.com

Distributors in other regions are listed on www.biotage.com

Literature Number: AN146

© 2020 Biotage. All rights reserved. No material may be reproduced or published without the written permission of Biotage. Information in this document is subject to change without notice and does not represent any commitment from Biotage. E&OE. A list of all trademarks owned by Biotage AB is available at www.biotage.com/legal. Other product and company names mentioned herein may be trademarks or registered trademarks and/or service marks of their respective owners, and are used only for explanation and to the owners' benefit, without intent to infringe.

