

ISOLUTE® cSPE for QuEChERS

ISOLUTE® cSPE for QuEChERS are pre-packed columns designed for efficient clean-up of pesticide residues in fruit and vegetables. The columns use the same sorbent blends as dispersive solid phase extraction (dSPE) QuEChERS, and are designed to replace the traditional dSPE step with column solid phase extraction (cSPE). The columns are available for both AOAC and EN methods and for general, waxed and pigmented fruits and vegetables.

Compared to traditional QuEChERS dSPE clean-up, these columns deliver:

- » A simplified workflow with fewer time-consuming manual steps
- » Cleaner extracts
- » Improved reproducibility
- » Automated processing options for improved accuracy

QuEChERS Clean-up Workflow

Typically, clean-up of pesticide residues in fruit and vegetables using QuEChERS methods relies on the dSPE workflow where the sample is shaken with dSPE media and then centrifuged to separate the cleaned-up sample. This is time consuming and complex because of the numerous manual transfer, capping and de-capping steps involved.



Figure 1. ISOLUTE® cSPE columns for QuEChERS

With pre-packed cSPE columns the workflow is significantly enhanced using a convenient flow-through format. So, instead of the usual procedure, the columns can be processed using vacuum, manual positive pressure or automated systems, leading to increased throughput. The simplified workflow means fewer manual transfer steps, which improves reproducibility, and reduces the risk of errors leading to lost samples or repeats.

The simplified cSPE process is illustrated in figure 2. Unlike dSPE, there is no need to shake the QuEChERS extract with dSPE media, or centrifuge to separate out the cleaned-up extract. Simply add the extract to the cSPE column, apply gentle vacuum or pressure, and collect the clean extract directly into a vial, ready for analysis.




Figure 2. Traditional dSPE clean-up workflow vs ISOLUTE® cSPE for QuEChERS clean-up workflow

Using a variety of produce types, including apples, oranges, oats, broccoli and spinach, we have compared the performance of traditional dSPE with the cSPE format. Analyte recovery of a representative suite of 56 pesticides is comparable between the two formats. Figure 5 shows pesticide recovery from oats using dSPE and cSPE clean-up. The enhanced clean-up using the column format means that some pesticides that are not quantifiable using dSPE can be successfully quantified using cSPE clean-up. In this case six pesticides were not detected using dSPE clean-up (myclobutanil, fenhexamide, tebuconazole, spiromesifen, bifentazate, fenoxycarb) whereas all were quantified in the cSPE extract.

Variability in analyte recovery from sample to sample is also improved in cSPE, due to the reduced number of manual steps, as well as the ease of quantitation due to cleaner extracts. This is demonstrated in figure 6 below. %RSD (n=7) across the pesticide panel for cSPE were significantly lower %RSDs than in the equivalent dSPE experiment.




Figure 6. cSPE provides superior reproducibility compared to dSPE. For a panel of 56 pesticides spiked into oranges, extracted using EN QuEChERS salts, and cleaned-up using ISOLUTE® EN Waxed media, 35 (~60%) had %RSD of <5%, compared to only 18 (~30%) with dSPE. In addition, five pesticides (captan 2, fenhexamide, tebuconazole, spiromesifen and acetamipride) were not detected using dSPE, with only one (spiromesifen) not detected using cSPE.

Automated cSPE Processing



Figure 7. Biotage® Extrahera™ HV-5000 automated sample preparation workstation

The ISOLUTE® cSPE format is ideal for automation of the QuEChERS clean-up step. Biotage® Extrahera™ HV-5000 automated sample preparation workstation (Fig. 7) can process a batch of 12 samples in just over 10 minutes. QuEChERS salt extracts are transferred from the homogenization/extraction step (performed in 50 mL tubes) directly into the sample rack without the need for any manual sample transfer.




Figure 8. Extrahera's intelligent pipetting can transfer the top QuEChERS layer cleanly from the QuEChERS homogenization/extraction tube directly to the cSPE clean-up column

Extrahera's intelligent pipetting feature means that QuEChERS extracts can be accurately loaded onto the cSPE column without perturbation of the sample or risking contamination of the clean-up column (Fig. 8). The entire volume of cleaned-up extract can then be collected directly into a vial suitable for subsequent evaporation or analytical steps.

Using smaller scale homogenization, or manual transfer from 50 mL homogenization tubes into sample tubes, batches of up to 48 samples can be cleaned-up in as little as 33 minutes*.

Standardized Media Blends

Standard AOAC and EN QuEChERS dSPE media blends are available pre-packed in this easy-to-use format, so you can match the clean-up to the appropriate sample matrix, delivering consistently effective clean-up and high analyte recoveries.

ISOLUTE® cSPE QuEChERS Media	QuEChERS Application
AOAC General	Clean-up of general fruit and vegetables (AOAC method) e.g. celery, head lettuce, cucumber, melon
AOAC Waxed	Clean-up of fruit and vegetables with wax or fat content (AOAC method) e.g. cereals, avocado seeds, nuts, dairy
AOAC Pigment	Clean-up of pigmented fruit and vegetables (AOAC method) e.g. strawberries, sweet potatoes, tomatoes
EN General	Clean-up of general fruit and vegetables (EN method) e.g. celery, head lettuce, cucumber, melon
EN Waxed	Clean-up of fruit and vegetables with wax or fat content (EN method) e.g. cereals, avocado seeds, nuts, dairy
EN Pigment	Clean-up of pigmented fruit and vegetables (EN method) e.g. strawberries, sweet potatoes, tomatoes
EN High Pigment	Clean-up of highly pigmented fruit and vegetables (EN method) e.g. red peppers, spinach, blueberries

* Using Biotage® Extrahera® Classic

Ordering Information

ISOLUTE® cSPE for QuEChERS Clean-up

Part Number	Description	Pack size
Q0030-0020-BG	ISOLUTE® AOAC General 200 mg/3 mL (Tablet)	50
Q0050-0035-BG	ISOLUTE® AOAC Waxed 350 mg/3 mL (Tablet)	50
Q0070-0022-BG	ISOLUTE® AOAC Pigment 225 mg/ 3 mL (Tablet)	50
Q0035-0020-BG	ISOLUTE® EN General 200 mg/3 mL (Tablet)	50
Q0060-0035-BG	ISOLUTE® EN Waxed 350 mg/3 mL (Tablet)	50
Q0080-0030-BG	ISOLUTE® EN Pigment 300 mg/3 mL (Tablet)	50
Q0090-0050-BG	ISOLUTE® EN High Pigment 500 mg/3 mL (Tablet)	50

EUROPE

Main Office: +46 18 565900
Toll Free: +800 18 565710
Fax: +46 18 591922
Order Tel: +46 18 565710
Order Fax: +46 18 565705
order@biotage.com
Support Tel: +46 18 56 59 11
Support Fax: + 46 18 56 57 11
eu-1-pointsupport@biotage.com

NORTH & LATIN AMERICA

Main Office: +1 704 654 4900
Toll Free: +1 800 446 4752
Fax: +1 704 654 4917
Order Tel: +1 704 654 4900
Order Fax: +1 434 296 8217
ordermailbox@biotage.com
Support Tel: +1 800 446 4752
Outside US: +1 704 654 4900
us-1-pointsupport@biotage.com

JAPAN

Tel: +81 3 5627 3123
Fax: +81 3 5627 3121
jp_order@biotage.com
jp-1-pointsupport@biotage.com

CHINA

Tel: +86 21 68162810
Fax: +86 21 68162829
cn_order@biotage.com
cn-1-pointsupport@biotage.com

KOREA

Tel: + 82 31 706 8500
Fax: + 82 31 706 8510
korea_info@biotage.com
kr-1-pointsupport@biotage.com

Distributors in other regions
are listed on www.biotage.com

Part Number: PPS719

© 2023 Biotage. All rights reserved. No material may be reproduced or published without the written permission of Biotage. Information in this document is subject to change without notice and does not represent any commitment from Biotage. E&OE. A list of all trademarks owned by Biotage AB is available at www.biotage.com/legal. Other product and company names mentioned herein may be trademarks or registered trademarks and/or service marks of their respective owners, and are used only for explanation and to the owners' benefit, without intent to infringe.