

# Maximize antibody recovery with PrismaA PhyTip® columns

• Gabrielle Kerkow, Application Scientist


## High capacity and broad range antibody purification

Biotage PrismaA PhyTip® columns leverage a unique combination of patented dual flow chromatography technology and Cytiva's MabSelect Prisma™ resin. This powerful pairing enables **high throughput, small volume, automated purification** while benefiting from superior binding capacity of the PrismaA resin. The key benefits of this approach are:

- **High performance:** Efficient and effective purification.
- **Minimum dilution:** Concentrated final product.
- **High throughput:** Process more samples in less time.
- **Small volume samples:** Ideal for limited sample sizes.
- **Scalable:** Easily adaptable for different volumes and needs.
- **Reproducible:** Consistent and reliable results.

## PrismaA PhyTip® Columns: Enhanced binding capacity with increased residence time

When purifying human IgG spiked in a 5x capture buffer, all three PhyTip® columns - ProA, ProPlus, and PrismaA - achieved similar recovery levels. However, PrismaA PhyTip® columns demonstrated a clear performance advantage, delivering a **30% higher yield** than ProPlus when the flow rate was reduced to 4 µL/sec (Figure 1).




**Figure 1.** Comparison of percent recovery of PhyTip® Columns. Reconstituted Human-IgG in 5x Capture Buffer at 0.2mg/mL starting concentration and 500 µL starting volume. Significant increase in yield is observed in PrismaA PhyTip® Column when flow rate is reduced.



## High-capacity serum antibody purification

Under identical purification conditions (an 8-cycle capture at a flow rate of 4  $\mu\text{L}/\text{sec}$ ), the ProPlus PhyTip® columns yielded 29% more antibodies from serum than the ProA columns. Impressively, the newer PrismaA PhyTip® columns demonstrated even higher capacity, delivering a **68% increase over ProA** and an additional **30% gain compared to ProPlus** (Figure 2).


PhyTip® column	Avg conc (mg/mL)	Avg yield ( $\mu\text{g}$ )	Capacity % increase from ProA	Capacity % increase from ProPlus
ProA	3.09	463.88	-	-
ProPlus	3.99	598.28	29%	-
PrismaA	5.19	778.85	68%	30%



**Figure 2.** Purification of antibodies from Porcine Serum using PhyTip® Columns. Serum was diluted 1:2 with 1x Capture Buffer. Starting volume was 1000  $\mu\text{L}$ . Samples were all processed with the same capture cycles, flow rate, and other conditions.



Further testing revealed that reducing the flow rate from 8  $\mu\text{L}/\text{sec}$  to 4  $\mu\text{L}/\text{sec}$  resulted in 20% higher yield with PrismaA PhyTip® columns (Figure 3). Notably, the same flowrate adjustment did not improve yields for ProPlus PhyTip® columns, underscoring PrismaA's superior binding efficiency and responsiveness to residence time.



**Figure 3.** Effects of Flow Rate on PhyTip® Columns. A significant increase in yield when reducing flow rate is observed in PrismaA PhyTip® Column from 8  $\mu\text{L}/\text{s}$  to 4  $\mu\text{L}/\text{s}$ . This is not observed in ProPlus PhyTip® Columns. Samples here were 1:2 Porcine Serum diluted with 1x Capture Buffer. Starting volume 500  $\mu\text{L}$ .

## Conclusion:

### Choosing between ProPlus and PrismaA: Speed vs. maximum capacity

If time is a limiting factor and a higher capacity than ProA is needed, ProPlus PhyTip® columns offer a strong performance advantage without requiring slower flow rates.

For users seeking maximum antibody recovery and binding capacity, PrismaA PhyTip® columns are the optimal choice. However, to unlock their full potential, a reduced flow rate during the capture step is essential. Notably, increasing the number of capture cycles did not significantly improve binding efficiency with PrismaA resin, highlighting that flow rate - not cycle count - is the key driver of performance.