

Biotage performance on the system you already own

The problem

- Prep HPLC/custom skids can be inconsistent at large scale.
- Replacing infrastructure is costly and slow.
- Manual packing limits reproducibility and throughput.

The Biotage answer

Plug and Play modules third-party integration — Bring FL150/FL400 radial compression modules in line with your existing pumps, hoists, and controls. Use prepacked, ready to run columns for stable flow and reproducible results, even in HP-API/cGMP workflows. Documentation (materials, pressure tests, schematics) included.



What you gain


- Rapid integration with minimal disruption.
- Higher throughput and better reproducibility vs. self-packed or high-pressure alternatives.
- Lower risk & reduced cycle times cross batches.

Best fit

Large scale HPLC or bespoke skids; CDMOs with validated infrastructure seeking Biotage level purification without major CAPEX.



Scale faster, run safer — on the system you already own.



Literature Number: PPS780.1

© 2026 Biotage. All rights reserved. No material may be reproduced or published without the written permission of Biotage. Information in this document is subject to change without notice and does not represent any commitment from Biotage. E&OE. A list of all trademarks owned by Biotage AB is available at www.biotage.com/legal. Other product and company names mentioned herein may be trademarks or registered trademarks and/or service marks of their respective owners, and are used only for explanation and to the owners' benefit, without intent to infringe.